




CHUNG INSTITUTE
OF INTEGRATIVE MEDICINE

DIRECTIONS FOR KMT

KMT Number: _____

KMT Duration Time: _____



1. Press and hold the ON button until you hear the beep followed by a flashing lights.
2. Rub alcohol on the leads prior to placing on the skin. Leads should be placed side by side by color: Red/Green – “christmas tree”, Yellow/Black – “bumble bee”.
3. After placing the proper leads on the assigned points, press SEL until you have selected the number you were given.
4. Once you have selected your number, press RUN button once. You should see the ON and KMT number flashing. This means the machine has started.